


Product Datasheet

CLCN6 antibody (orb543719)

Description	Rabbit polyclonal antibody to CLCN6
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human CLCN6(NP_001277.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51797
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue showing CLCN6 expression.