



Product Datasheet

Doublecortin antibody (orb543730)


Description	Rabbit Polyclonal antibody to Doublecortin
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Doublecortin(NP_835365.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41KDa
Uniprot ID	O43602
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC-12 cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of HeLa cell...