



Product Datasheet

DNMT1 antibody (orb543735)


Description	Rabbit polyclonal antibody to DNMT1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human DNMT1(NP_001370.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200KDa
Uniprot ID	P26358
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse spinal cord tissue.



Immunohistochemical staining of human placental tissue.



Immunohistochemical staining of human colon tissue.