



Product Datasheet

PIK3R2 antibody (orb543778)


Description	Rabbit polyclonal antibody to PIK3R2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PI3 Kinase p85 beta (NP_005018.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	000459
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse spleen...



Immunohistochemical staining of mouse kidney...



Immunohistochemical staining of rat spleen...