



Product Datasheet

MonoMethyl-Histone H2B-K5 antibody (orb543811)


Description	Rabbit polyclonal antibody to MonoMethyl-Histone H2B-K5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic monomethylated peptide around K5 of human H2BK5me1(NP_003519.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	14kDa
Uniprot ID	Q16778
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunofluorescence analysis of U-2 OS ce...



Immunofluorescence analysis of NIH-3T3 C...