


## Product Datasheet

### TOMM7 antibody (orb543960)

|                            |                                                  |
|----------------------------|--------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TOMM7              |
| <b>Reactivity</b>          | Mouse                                            |
| <b>Conjugation</b>         | Unconjugated                                     |
| <b>Tested Applications</b> | WB                                               |
| <b>Immunogen</b>           | A synthetic peptide of human TOMM7(NP_061932.1). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.      |
| <b>Note</b>                | For research use only                            |
| <b>Isotype</b>             | IgG                                              |
| <b>Clonality</b>           | Polyclonal                                       |
| <b>MW</b>                  | 11kDa                                            |
| <b>Uniprot ID</b>          | <b>Q9P0U1</b>                                    |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |



Western blot analysis of extracts of Mou...