


Product Datasheet

TUBGCP5 antibody (orb544061)

Description	Rabbit polyclonal antibody to TUBGCP5
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human TUBGCP5(NP_443135.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	118kDa
Uniprot ID	Q96RT8
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Rat...