



Product Datasheet

BMPR1B antibody (orb544103)


Description	Rabbit polyclonal antibody to BMPR1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant protein of human BMPR1B.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	60kDa
Uniprot ID	O00238
Entrez	658
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human ma...