


Product Datasheet

MonoMethyl-Histone H4-K5 antibody (orb544104)

Description	Rabbit polyclonal antibody to MonoMethyl-Histone H4-K5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	A synthetic peptide of human MonoMethyl-Histone H4-K5.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P62805
Entrez	8370
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...