

## Product Datasheet

### PRKCA antibody (orb544124)

|                            |                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PRKCA                                                                                                           |
| <b>Reactivity</b>          | Mouse                                                                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                  |
| <b>Tested Applications</b> | WB                                                                                                                                            |
| <b>Immunogen</b>           | A synthetic peptide of human Pan PKCPKC.                                                                                                      |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.                                                                                               |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                                                    |
| <b>MW</b>                  | 83kDa                                                                                                                                         |
| <b>Uniprot ID</b>          | <b>Q02156</b> , <b>P17252</b> , <b>P05771</b> , <b>Q05513</b> , <b>P41743</b> , <b>Q05655</b> , <b>Q04759</b> , <b>P05129</b> , <b>P24723</b> |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                               |