



Product Datasheet

M6A antibody (orb544126)

Description	Rabbit polyclonal antibody to M6A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Chemical compounds corresponding to N6-methyladenosine / m6A.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Dot-blot analysis of Oligomer 1 - N6-met...



Dot-blot analysis of M6A antibody are te...