

Product Datasheet

ETV5 antibody (orb544129)

Description	Rabbit polyclonal antibody to ETV5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ETV5(P41161).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62kDa
Uniprot ID	P41161
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.