

Product Datasheet

ADAMTS9 antibody (orb544130)

Description	Rabbit polyclonal antibody to ADAMTS9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human ADAMTS9.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9P2N4
Entrez	56999
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.