

Product Datasheet

CTHRC1 antibody (orb544168)

Description	Rabbit polyclonal antibody to CTHRC1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CTHRC1(NP_612464.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24kDa/25kDa/26kDa
Uniprot ID	Q96CG8
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.