

Product Datasheet

SLC30A7 antibody (orb544171)

Description	Rabbit polyclonal antibody to SLC30A7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLC30A7 (NP_598003.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q8NEW0
Entrez	148867
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.