


Product Datasheet

CNOT1 antibody (orb544179)

Description	Rabbit polyclonal antibody to CNOT1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CNOT1(NP_057368.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	267kDa
Uniprot ID	A5YKK6
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...