

Product Datasheet

Galectin 3 antibody (orb544196)

Description	Rabbit Polyclonal antibody to Galectin 3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human LGALS2(NP_006489.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15KDa
Uniprot ID	P05162
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.