



Product Datasheet

Acetyl-Histone H4-K91 antibody (orb544214)


Description	Rabbit polyclonal antibody to Acetyl-Histone H4-K91
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic Peptide of human Acetyl-Histone H4-K91
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	17kDa
Uniprot ID	P62805
Entrez	8370
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:200, ChIP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of C6 cells ...