

Product Datasheet

NLE1 antibody (orb544221)

Description	Rabbit polyclonal antibody to NLE1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NLE1(NP_060566.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa/53kDa
Uniprot ID	Q9NVX2
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.