

Product Datasheet

KLHL7 antibody (orb544235)

Description	Rabbit polyclonal antibody to KLHL7
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human KLHL7(NP_061334.4).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q8IXQ5
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.