

Product Datasheet

P2RX1 antibody (orb544241)

Description	Rabbit polyclonal antibody to P2RX1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human P2RX1 (NP_002549.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P51575
Entrez	5023
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.