

Product Datasheet

SEPT6 antibody (orb544261)

Description	Rabbit polyclonal antibody to SEPT6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Septin 6(NP_055944.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa
Uniprot ID	Q14141
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.