

Product Datasheet

TEKT4 antibody (orb544304)

Description	Rabbit polyclonal antibody to TEKT4
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TEKT4(NP_653306.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q8WW24
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.