

Product Datasheet

EYA1 antibody (orb544308)

Description	Rabbit polyclonal antibody to EYA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human EYA1 (NP_742055.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	55kDa
Uniprot ID	Q99502
Entrez	2138
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.