

Product Datasheet

Human ISCU protein (orb54431)

Description	Recombinant human ISCU protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	N-terminal GST-tagged: N-terminal GST-tagged1-308AA: 1-167AFull Length : Full Length
Protein Sequence	YHKKVVDHYENPRNVGSLDKTSKNVGTGLVGAPACGDVMKLQIQVDEKGIKVD ARFKTFGCGSAIASSSLATEWVKGKTVEEALTIKNTDIAKELCLPPVKLHCSMLAE DAIKAALADYKLKQEPKKGEAEKK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	18.0 kDa
Uniprot ID	Q9H1K1
Expiration Date	6 months from date of receipt.