

Product Datasheet

GNGT2 antibody (orb544318)

Description	Rabbit polyclonal antibody to GNGT2
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GNGT2(NP_113686.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	O14610
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.