

Product Datasheet

KCNJ9 antibody (orb544319)

Description	Rabbit polyclonal antibody to KCNJ9
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human KCNJ9(NP_004974.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	Q92806
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.