

Product Datasheet

INSR (Phospho-Y1146)+ IGF1R (Phospho-Y1131) antibody (orb544323)

Description	Rabbit polyclonal antibody to INSR (Phospho-Y1146)+ IGF1R (Phospho-Y1131)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1146 of human INSR and Y1131 of human IGF1R
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P08069
Entrez	34803643
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.