

Product Datasheet

MAP2K3 (Phospho-S189)+ MAP2K6 (Phospho-S207) antibody (orb544328)

Description	Rabbit polyclonal antibody to MAP2K3 (Phospho-S189)+ MAP2K6 (Phospho-S207)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S189 of human MAP2K3 and S207 of human MAP2K6
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	40KDa
Uniprot ID	P46734
Entrez	56065608
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.