


Product Datasheet

RB1 (Phospho-S780) antibody (orb544331)

Description	Rabbit polyclonal antibody to RB1 (Phospho-S780)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S780 of human RB
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	100kDa
Uniprot ID	P06400
Entrez	5925
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NIH...