

Product Datasheet

BRCA1 (Phospho-S1524) antibody (orb544341)

Description	Rabbit polyclonal antibody to BRCA1 (Phospho-S1524)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S1524 of human BRCA1 (NP_009225.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P38398
Entrez	672
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.