

Product Datasheet

CHUK (Phospho-S176) antibody (orb544342)

Description	Rabbit polyclonal antibody to CHUK (Phospho-S176)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S176 of human IKK α (NP_001269.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	84kDa
Uniprot ID	O15111
Entrez	1147
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.