

Product Datasheet

ESR1 (Phospho-S118) antibody (orb544343)

Description	Rabbit polyclonal antibody to ESR1 (Phospho-S118)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S118 of human Estrogen Receptor alpha (ESR1) (NP_000116.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P03372
Entrez	2099
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.