

Product Datasheet

NFKB1 (Phospho-S933) antibody (orb544346)

Description	Rabbit polyclonal antibody to NFKB1 (Phospho-S933)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S933 of human NFKB1 (NP_001158884.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P19838
Entrez	4790
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.