

Product Datasheet

PLK1 (Phospho-T210) antibody (orb544349)

Description	Rabbit polyclonal antibody to PLK1 (Phospho-T210)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T210 of human PLK1(NP_005021.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	68kDa
Uniprot ID	P53350
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.