

Product Datasheet

CSF1R (Phospho-Y699) antibody (orb544373)

Description	Rabbit polyclonal antibody to CSF1R (Phospho-Y699)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y699 of human CSF1R (NP_005202.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P07333
Entrez	1436
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.