

Product Datasheet

CSF1R (Phospho-Y708) antibody (orb544375)

Description	Rabbit polyclonal antibody to CSF1R (Phospho-Y708)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y708 of human CSF1R(NP_005202.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	175KDa
Uniprot ID	P07333
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.