

Product Datasheet

Human BCL2L15 protein (orb54438)

Description	Recombinant human BCL2L15 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	N-terminal GST-tagged: N-terminal GST-tagged1-206AA: 1-163AAFull Length : Full Length of BC127719
Protein Sequence	MKSSQTFEEQTECIVNTLLMDFLSPTLQVASRNLCCVDEVDSGEPCSFVDVAIIAG RLRMLGDQFNGELEASAKNVIAETIKGTGAILQNTVESLSKTWCAQDSSLAYER AFLAVSVKLLEYMAHIAPEVVGQVAIPMTGMINGNQAIREFIQGGGWENLES
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	44.7 kDa
Uniprot ID	Q5TBC7
Expiration Date	6 months from date of receipt.