

## Product Datasheet

### ARAF (Phospho-S299) antibody (orb544380)

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ARAF (Phospho-S299)                          |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                          |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | WB                                                                         |
| <b>Immunogen</b>           | A synthetic phosphorylated peptide around S299 of human ARAF(NP_001645.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.                           |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                |
| <b>Note</b>                | For research use only                                                      |
| <b>Isotype</b>             | IgG                                                                        |
| <b>MW</b>                  | 20kDa/67kDa                                                                |
| <b>Uniprot ID</b>          | <b>P10398</b>                                                              |
| <b>Dilution Range</b>      | WB: 1:1000-1:2000                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |