


Product Datasheet

ULK1 (Phospho-S757) antibody (orb544384)

Description	Rabbit polyclonal antibody to ULK1 (Phospho-S757)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S757 of human Phospho-ULK1-S757 (NP_003556.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	145KDa
Uniprot ID	075385
Entrez	8408
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A431...