

Product Datasheet

CAV1 (Phospho-Y14) antibody (orb544385)

Description	Rabbit polyclonal antibody to CAV1 (Phospho-Y14)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y14 of human Caveolin-1 (NP_001744.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	23kDa
Uniprot ID	Q03135
Entrez	857
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.