

Product Datasheet

eIF4E (Phospho-S209) antibody (orb544386)

Description	Rabbit polyclonal antibody to eIF4E (Phospho-S209)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S209 of human eIF4E(NP_001959.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	25kDa
Uniprot ID	P06730
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.