

Product Datasheet

ULK1 (Phospho-S555) antibody (orb544387)

Description	Rabbit polyclonal antibody to ULK1 (Phospho-S555)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S555 of human Phospho-ULK1-S555 (NP_003556.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	113KDa
Uniprot ID	O75385
Entrez	8408
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.