

Product Datasheet

PXN (Phospho-Y118) antibody (orb544389)

Description	Rabbit polyclonal antibody to PXN (Phospho-Y118)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y118 of human Paxillin (NP_001074324.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	56kDa
Uniprot ID	P49023
Entrez	5829
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.