

Product Datasheet

GAB1 (Phospho-Y627) antibody (orb544392)

Description	Rabbit polyclonal antibody to GAB1 (Phospho-Y627)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y627 of human GAB1 (NP_002030.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q13480
Entrez	2549
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.