

Product Datasheet

eIF4B (Phospho-S406) antibody (orb544395)

Description	Rabbit polyclonal antibody to eIF4B (Phospho-S406)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S406 of human eIF4B (NP_001408.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P23588
Entrez	1975
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.