

Product Datasheet

MAP3K7 (Phospho-T187) antibody (orb544400)

Description	Rabbit polyclonal antibody to MAP3K7 (Phospho-T187)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T187 of human TAK1(NP_003179.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	53kDa/56kDa/64kDa/67kDa
Uniprot ID	O43318
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.