

Product Datasheet

FOXO3 (Phospho-S294) antibody (orb544412)

Description	Rabbit polyclonal antibody to FOXO3 (Phospho-S294)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S294 of human FOXO3A(NP_001446.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	82kDa
Uniprot ID	O43524
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.