

Product Datasheet

STAT1 (Phospho-Y701) antibody (orb544414)

Description	Rabbit polyclonal antibody to STAT1 (Phospho-Y701)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y701 of human STAT1.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	83kDa
Uniprot ID	P42224
Entrez	6772
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.