

Product Datasheet

TP53 (Phospho-S392) antibody (orb544416)

Description	Rabbit polyclonal antibody to TP53 (Phospho-S392)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S392 of human p53(NP_000537.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	53kDa
Uniprot ID	P04637
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.